CANON MEDICAL SYSTEMS DO BRASIL LTDA.

A Ultrassonografia Pulmonar na Covid-19

Aline de Gois Collete

Dr. Thiago Potrich

Dra. Maria Cristina Chammas

Introdução

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
Score	Findings
Score 0 – normal aeration	A lines – maximum 2 B-lines
Score 1 – Moderate loss of aeration	Artifacts occupying ≤ 50% of the pleura
Score 2 – Severe loss of aeration	Artifacts occupying > 50% of the pleura
Score 3 – Complete loss of aeration	Tissue-like pattern (lung consolidation)

Figura 12 *Lung Ultrasound Score – LUS* - de Alencar, J.C.G., Marchini, J.F.M., Marino, L.O. et al. Lung ultrasound score predicts outcomes in COVID-19 patients admitted to the emergency department. *Ann. Intensive Care* **11**, 6 (2021)

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1. Oliveira RR, Rodrigues TP, Savoia P, Gomes AC, Chammas MC. Lung Ultrasound: na additional tool in COVID-19. *Radiol Bras.* 53 (4) 241-251, 2020.
2. Lichtenstein DA, Mezière GA. Relevance of lung ultrasound in the diagnosis of acute respiratory failure: the BLUE protocol. *Chest.* 2008;134:117-25.
3. Buonsenso D, Pata D, Chiaretti A. COVID-19 outbreak: less stethoscope, more ultrasound. *Lancet Respir Med.* 2020;8:e27.
4. Soldati G, Smargiassi A, Inchingolo R, et al. Proposal for international standardization of the use of lung ultrasound for patients with COVID-19: a simple, quantitative, reproducible method. *J Ultrasound Med.* 2020 Mar 30.
5. Verbeeh JH, Rajamaki B, Ijaz S, Sauni R, Toomey E, Blackwood , Tikka C, et al. Personal protective equipment for preventing highly infectious diseases due to exposure to contaminated body fluids in healthcare staff. *Cochrane Database of Systematic Reviews.* Issue 5. Art. No.: CD011621, 2020.
6. Koh DM, Burke S, Davies N, Padley S P. Transthoracic US of the Chest: Clinical Uses and Applications. *RadioGraphics* 22 (1), 2002.
7. Koh DM, Burke S, Davies N, Padley S P. Transthoracic US of the Chest: Clinical Uses and Applications. *RadioGraphics* 22 (1), 2002.
8. Farias LP, Fonseca EK, Strabelli DG, Loureiro BM, Neves YC, Rodrigues TP, et al. Imaging findings in COVID-19 pneumonia. *Clinicis.* São Paulo. *Clinicis.* Vol 75:e2027, 2020.
9. Clevert DA. White paper: Lung Ultrasound in Patients with Coronavirus COVID-19 Disease. [siemens-healthineers.com/ultrasound](https://www.siemens-healthineers.com/ultrasound). 2020.
10. Ng MY, Lee EYP, Yang J, et al. Imaging profile of the COVID-19 infection: radiologic findings and literature review. *Radiology: Cardiothoracic Imaging.* 2020;2(1):e200034.
11. Gargani L, Volpicelli, G. How I do it: Lung ultrasound. *Cardiovasc Ultrasound.* 12 (25), 2014.
12. Lichtenstein DA, Mezière GA, Lagoueyte JF, et al. A-lines and B-lines: lung ultrasound as a bedside tool for predicting pulmonary artery occlusion pressure in the critically ill. *Chest.* 2009;136:1014-20.
13. Lichtenstein D, Mezière G, Seitz J. The dynamic air bronchogram. A lung ultrasound sign of alveolar consolidation ruling out atelectasis. *Chest.* 2009;135:1421-5.

Canon
CANON MEDICAL SYSTEMS DO BRASIL

br.medical.canon

©Canon Medical Systems Brasil 2020. All rights reserved.

UL Xario 200G

Made for Life are trademarks of Canon Medical Systems Corporation

Made For life